

Rapid Test for Detection of Antigen of SARS-CoV-2

CareStart™ COVID-19 Antigen Home Test

 **US FDA**
Emergency
Use Authorized



The CareStart™ COVID-19 Antigen Home Test is a lateral flow immunoassay intended for the qualitative detection of nucleocapsid protein antigens from SARS-CoV-2.

This test is authorized for non-prescription home use with self-collected anterior nasal (nares) swab samples from individuals aged 14 years or older or adult collected anterior nasal (nares) swab samples from individuals aged 2 years or older.

- This test is authorized for individuals with symptoms of COVID-19 within the first 7 days of symptom onset when tested at least twice over three days with at least 48 hours between tests, and for individuals without symptoms or other epidemiological reasons to suspect COVID-19, when tested at least three times over five days with at least 48 hours between tests.
- For use under Emergency Use Authorization (EUA) only
- For *in vitro* diagnostic use only

Home Test for SARS-CoV-2 Antigen Detection (Ages 2 and Up)

- ✓ Fast and easy to self-test anywhere
- ✓ Easy to use with mobile application providing step-by-step instructions
- ✓ Test results obtained in 10 minutes
- ✓ Identify individual's current infection status to COVID-19
- ✓ No special equipment or training required
- ✓ Health Canada authorized for use under Interim Order