


M-Series Short Infection Control Cart,

Mini Width, Four Drawers, E-Lock,

M3DS1824E03



Available Colors



Harloff's M3DS1824E03 offers customers an ideal infection control cart that can easily be stored in any facility, office or room where space is limited. Constructed from heavy-duty steel and featuring an electronic pushbutton lock with access for up to 20 users, this infection control cart can be customized with multiple locking options depending on your needed level of security. This infection control cart also features 24" of vertical drawer space with customizable drawer configurations, full extension drawer slides, a removeable antimicrobial plastic top and 3" single wheel stem casters with two braking casters. The M3DS1824E03 replaces the Classic model 3143E.

Features

- Sturdy, uni-body steel construction
- Basic electronic pushbutton lock
- Key lock override
- Factory set user and supervisor code
- Cabinet may be programmed in the field
- Durable scratch resistant heat sealed powder coat finish
- Full extension drawers slides for easy access
- Body has pre-threaded accessory mounting holes
- Antimicrobial plastic top, removable
- 24" of vertical drawer space
- Cart is delivered assembled
- Fully customizable drawer configurations
- 3" (7.6 cm) single wheel stem caster set, full swivel - two with brake
- Latex free
- 12-Year Manufacturer's Warranty on shipments to US and Canada. 1-Year Limited Warranty on any bought-out parts.

Email us:
sraithel@products-unlimited.com

Visit us:
products-unlimited.com

Contact us:
800-865-4683