

# M3 Vital Sign Monitor

Part Number: ED-M3NSTO-EA



**User-friendly. Sleek. The M3 serves as a vital spot check and a continuous vital signs monitor providing a comprehensive solution for medical professionals.**

## HIGHLIGHTED FEATURES:

- 5.6" High resolution color TFT-LCD screen
- Lightweight & portable design
- Spot and continuous monitoring
- Flexible configurations to meet different clinical needs
- Auto/Manual/Continuous/Average BP mode
- Real-time measurements & trend display
- Built-in rechargeable Lithium-ion battery for 8 working hours
- User-friendly interface
- Nurse call
- Optional thermal printer



## COMMUNICATION FEATURES:

- Bi-directional communication with MFM-CMS central monitoring system
- Compatible with CMS Lite PC data management software
- HL7/AIP
- Wired LAN/WiFi Support
- USB port for flash drive or pocket router for optional WiFi
- Barcode scanner support
- Export trends as XML or PDF in the central monitoring system.