

Pink Adscope Stethoscope

SKU: AD-603PBCA-EA



DESCRIPTION:

One of our most popular acoustic scopes, the ADC Adscope 603 combines outstanding performance with rugged durability. The combination chestpiece offers both an ultra-sensitive diaphragm for greater amplification and crisper high-frequency transmission as well as an extra-deep bell for unsurpassed low-frequency response. Our adjustable frequency (AFD) technology allows you to tune the diaphragm to further enhance high- and low-frequency response. The clinician headset, designed to withstand thousands of cycles, features a double-leaf tensioning spring and a reinforcing yoke molded into flexible 22-inch PVC tubing. The Adscope 603 features an all stainless-steel construction, ensuring a lifetime of use (tactical finish has black anodized aluminum binaurals).